

SMARTeZ® Infusionspumpen

kurzfristige Infusionen



Tabelle der Auswirkung von einer grösseren oder geringeren Füllmenge

Model Nr (REF)	481030	481090	481150	481040	481120	481100	481010	481050	481080	481020	481060
Durchfluss (ml/st)	50	50	50	100	100	100	200	250	250	500	175
Auffüllungsvolumen (ml)	50	100	250	100	250	400	100	250	500	250	250
Minimal Volumen (ml)	50	75	150	50	200	275	50	150	360	250	250
Maximal Volumen(ml)	60	110	275	120	300	550	100	300	550	300	300
Restvolumen (ml)	2	2	3	2	3	5	2	3	5	3	3

Dauer der Infusion											
Stunden	Auffüllungsvolumen (ml)										
00:15							50(*)				
00:30				50			100	150		250	
00:45				75				200		300	140
01:00	50			100				250			175
01:15	60			120				300	360		220
01:30		75							400		250
01:45									475		300
02:00		100			200				500		
02:15		110			220				550		
02:30					250						
02:45					275	275					
03:00			150		300	300					
03:30						360					
04:00						400					
04:30						450					
05:00			250			500					
05:30			275			550					