

Tabelle der Auswirkung von einer grösseren oder geringeren Füllmenge

Model Nr ( REF)	480050	480080	480010	480030	480250	480090	480260	481200	481170	481180
Durchfluss (ml/st )	2	2	5	5	10	5	10	20	20	40
Auffüllungsvolumen (ml)	60	100	60	125	120	270	240	250	500	500
Minimal Volumen (ml)	50	75	60	90	90	250	200	200	350	240
Maximal Volumen(ml)	60	120	100	175	170	360	300	300	500	500
Restvolumen (ml)	2	2	2	2	480250	3			5	5

Dauer der Infusion										
Stunden	Tagen	Auffüllungsvolumen (ml)								
6 St										240
8 St						80				320
10 St						100		200		400
12 St				60		120		250	250	500
15 St								300		
18 St				90	90				360	
20 St				100			200		420	
25 St	1 T	50			125		240		500	
27 St										
30 St		60					300			
36 St	1,5 T		75		175					
50 St	2 T		100				250			
54 St			108				270			
72 St	3 T						360			